

STERILIZATION PORTAL

USER MANUAL



Company Information

Sterismart

4181 Sladeview Crescent Unit #36, Mississauga, L5L 5R2

Ontario, Canada

Office: (905) 997-7171

Email: info@sterismart.ca

Website: www.sterismart.ca

User Manual Version 1.00

Table of Contents

1	Introduction	5
1.1 <i>Sterismart</i>	5
1.1.1	Documented safety	5
1.2 <i>Sterilization Portal</i>	5
2	Getting Started	7
2.1 <i>System Setup</i>	7
2.2 <i>Registration</i>	7
2.3 <i>Registration Login</i>	7
2.4 <i>First Login</i>	10
3	Administration.....	11
3.1 <i>Administration Login</i>	11
3.2 <i>Setup/Change Office Information</i>	12
3.3 <i>View Support License Information</i>	13
3.4 <i>Setup Autoclave Information</i>	14
3.4.1	Add Autoclave	14
3.4.2	Missing Autoclave.....	16
3.4.3	Edit Autoclave	16
3.4.4	Delete Autoclave	16
3.5 <i>Setup Operator Information</i>	17
3.5.1	Add Operator	17
3.5.2	Edit Operator	18
3.5.3	Delete Operator.....	19
3.6 <i>Manage Packages</i>	20
3.6.1	Creating New Package	20
3.6.2	Modifying Existing Packages	21
3.6.3	Viewing Packages	22
3.7 <i>Backup and Restore</i>	22

3.7.1	Creating a Backup	22
3.7.2	Setting Up Automatic Backup	22
3.7.3	Restoring Data.....	23
3.8 Audit	24
3.9 Logging out from Administration Panel	27
4	Main Menu	28
4.1 Login.....	28
4.2 Navigation.....	28
4.2.1	Step 1 – Printing of Label	29
4.2.2	Step 2 – Sterilization.....	30
4.2.2.1	Smart Reader Data Logger	31
4.2.3	Step 3 – Biological Indicator Test	33
4.2.4	Step 4 – Reports	34
4.2.5	Administration.....	38
5	Help Center	39
5.1 Help.....	39
5.2 Training	39
5.3 FAQ.....	39
5.4 Documentation	39
5.5 Manual.....	39
6	Sign Out	39
7	Appendix	40
7.1 License and Serial Number	40

1 Introduction

1.1 Sterismart

Sterismart is part of the Hatch group of companies and is the software department. The sterilization portal is part of the Sterilization Suite built by Sterismart. For your patients, your team and your clinic!

The sterilization portal is one stop for all sterilization logs. Cycle number, temperature and pressure, lot number and type of load, biological indicator test and additional data in one software. Log cycle data and create unique barcode for every cycle. Scan the barcode and attach to client file in your patient management system.

1.1.1 Documented safety

The infection prevention and control, compliance and new technologies have changed the landscape of today's dental procedures. We have created the tool that will enhance and enable resources to complete the sterilization process under compliant rules. The Sterilization portal is the first choice for safe and professional documentation across the whole instrument treatment process. All process steps, logs and data are fully documented and archived securely. All logs and data entries are fully documented and archived securely. A report is created and sent as a tamper-proof PDF.

This will enable the dental office to ensure all logs and data are following regulations by public health, stored securely and provided by request to the public health authorities and patient. Barcode labels can then be printed and placed on the loads to facilitate traceability of the individual load elements. This can be read-in via a scanner and added to the patient administration system.

1.2 Sterilization Portal

The Sterilization Portal helps simplify sterilization data logging and management of log data as it compiles daily logs on autoclave cycles and biological indicator (B. I.) tests.

All in one place. Secure. In the local storage and designed to increase the efficiency and effectiveness of the workflows in dental offices of all sizes.

Sterismart delivers more functionality and value than our competition by the ability to use data from all autoclaves used by healthcare professionals.

- **Secure**

The login is secured by office ID and PIN. All reports contain all operators involved in the cycle.

- **Cycle Data, B. I. Test Report**

On completion of the cycle, you enter all the necessary data in the cycle log. Parameters are cleaning method, package ID(s), packages (name, type and number), autoclave name, cycle

number, temperature, pressure and chemical test results. The Sterilization Portal will store your results and attach them to the cycle report for the autoclave. The B. I. report data is stored and attached to the cycle logs of the autoclave. The Bowie Dick test for class B autoclaves is also recorded.

- **Barcode Tracking**

The system generates a barcode (or you can use pre-printed barcode labels) and your entry is complete. Simple, easy and quick. Print the barcode labels and apply to your sterilized packaging and store if not used immediately. Add the barcode to the patient file.

- **Reporting**

Sterismart creates sterilization reports that are compliant with Health Regulations and provides traceability report of sterilization cycle records. The Sterilization Portal will store your results and attach them to the cycle report for the autoclave.

The options for the reports are:

Cycle Report. You can generate cycle reports daily, by date Range, by autoclave, by sterilization user, by biological indicator test user or by package ID (barcode).

All reports are formatted to simplify the overview and have initials and signatures of operators and Managers.

2 Getting Started

2.1 System Setup


1. Please place the tablet in the position and place you wish to use
2. Please have a power outlet close to the unit
3. Please power on the system
4. Please click on the icon for the Sterismart software.

2.2 Registration

To start using the Sterilization Portal, you need to register your office. Please go to: <https://sterilizationportal.sterismart.ca>. You will be taken to the sterilization portal where you can complete your registration

2.4 Registration Login

1. From the documentation provided, you will find the serial number and license number to be used for your office set up.
2. Enter your serial number and license number.
3. Once you enter the serial number and the license number click on the continue button to login.



The screenshot shows the 'STERILIZATION PORTAL Setup Portal' interface. It features a teal background with a geometric pattern. At the top left, it says 'STERILIZATION PORTAL Setup Portal'. At the top right is the 'STERISmart' logo with the tagline 'Intelligent Log Management Solutions'. The main content area has the heading 'Please enter the information below to continue to set up your account'. Below this are two input fields: 'Please enter serial number' with a placeholder 'Enter serial number here' and 'Please enter license number' with a placeholder 'XXXX-XXXX-XXXX-XXXX'. A red arrow points from the text 'Serial number' to the first input field, and another red arrow points from the text 'License number' to the second input field. Below the input fields is a blue 'Continue' button. In the bottom right corner, there is a small blue square icon with a white speech bubble.

4. After successfully logging in, you will be taken to the office setup page. Please enter all the information requested.

5. Enter office name, example Smart Dental Mississauga.

STERILIZATION PORTAL

Office Setup

Please enter information below to set up office account

All * fields are required

Office Name *

Enter office PIN * Please enter 4 digits PIN of your choosing. This will be the office PIN for the office.

Administrator's First Name*

Administrator's Last Name*

Create Admin PIN * Please enter 4 digits PIN of your choosing. This will be the official PIN for the Admin.

Street Address * Unit Number

City * Postal Code *

Province * Country *

Enter Office Email * Telephone *

[Back to login](#) [Submit](#)

6. Enter a 4-digit pin to be used as office pin for login purpose. This pin will be shared with operators. Example 1234. Please make it a number that is easily remembered.
7. Please enter the Administrator First Name and Last Name.
8. Create the admin pin. This pin is NOT to be shared with operators. This is the admin pin. It will be used to access the administrative functions of Sterismart.
9. Please enter the street address, the unit number, city and the postal code, Province and the Country.
10. Please enter the main office email and telephone number for communication.
11. When all data is entered, please review and if all good, click on Submit.
12. Upon successful submission, you will get your office ID. This will send an email with information submitted in the form.

STERILIZATION PORTAL

Office Setup

Please enter information below to set up office account

All * fields are required

Office Name *

Enter office PIN * Please enter 4 digits PIN of your choosing. This will be the office PIN for the office.

Administrator's First Name*

Administrator's Last Name*

Create Admin PIN * Please enter 4 digits PIN of your choosing. This will be the official PIN for the Admin.

Street Address * Unit Number

City * Postal Code *

Province * Country *

Enter Office Email * Telephone *

[Back to login](#) [Submitting...](#)


Office account created

Your office ID is: LCHL1900

Please use this ID to log in to the application. You will receive an email with all log in credentials. Please proceed to initialize your desktop account.

[Continue](#)

13. Please note the office ID, as this is the primary login ID you will need. This information is also in the email.
14. You will then returned to the Setup Portal login page.
15. You will receive the confirmation email from the system.



The confirmation email will send the following information:

Hello Agobin Fobella,

Thank you for setting up your Sterismart account. Please find below your login information.

Office ID: LCHL1900,
Office PIN: 1234,
Administrator: Agobin Fobella,
Administrator's PIN: 1234.


Your account is now ready for use, please proceed to initialize your desktop account.
Sterismart Team.

Please print the email and consider it confidential. You will need the Office ID and Office PIN to log into the application. The administrator's PIN is for the administrator for setup and changes, as well as for work with the application.

2.5 First Login

Before you login in the desktop application, you will need to initialize your office. To do this, enter your office ID and PIN, you will be prompted, if you want to initialize your office, select yes. Office initialization will start. Upon successful initialization, you will be asked to login again to complete set up.

1. Please use your office ID and office PIN number to log into the application. Enter your office ID and PIN in their respective boxes and click on Login.




2. Upon successful login, you will be redirected to the administrator setup of the application where you will configure the software based on your office requirements.

3 Administration

3.1 Administration Login

The administrator can log into the admin login to perform changes, updates and deletes.

1. Click on the Administration button
2. Popup with request for PIN will show
3. Enter your administrator PIN
4. After a successful login, you will be redirected to the administration section



3.2 Setup/Change Office Information

In this section, you can modify your office information.

1. Please review office information
2. Click on the field you wish to change
3. After making changes, click on the Update Information button to update your office information.
4. Upon successful update, you will receive a notification that all changes have been successfully saved.

The screenshot displays the 'STERILIZATION PORTAL Dental One' web application. The top navigation bar includes a 'SIGN OUT' button, the 'STERI SMART' logo, and the date/time 'Thu, 09 Jan 2020 11:28'. A left sidebar contains a menu with options: LOGOUT, HELP CENTER, ABOUT US, OFFICE SETUP (highlighted), AUDIT, LICENSE INFO, AUTOCLAVE SETUP, OPERATOR SETUP, MANAGE PACKAGES, and BACKUP & RESTORE. The main content area is titled 'OFFICE SETUP' and contains a form with the following fields:

Office Name	Dental One		
Office ID	LCHL1900	Office PIN	1234
License Number	XXXX-XXXX-XXXX-0000	SN of Device	0000000000
Street Address	Street One	Unit Number	3
City	Brampton	Postal Code	L7U4I9
Province	Ontario	Country	Canada
Office Email	agobin@smartreader.ca		
Telephone	905-997-7171		
Report Email	<input type="text"/>		

Below the form is an 'Update Information' button. A red message at the bottom of the form states: 'Please enter report email to receive pdf reports'.

3.3 View Support License Information

With a valid support license, you will be able to get support for any issues that arise when you use the software. When your license expires, you will be required to renew it to continue receiving support.

1. This page shows the support license information of your office
2. You cannot change any information on this page
3. Once your license is 30 days or less from expiring, you will receive a notification and will be reminded by the app for renewal.
4. When your license expires, you will be directed to the renew license page where you will be able to renew your license.
5. You can also renew your support license from within the application. When your license is 30 day or less from expiring, the **“Click here to renew support license”** becomes active.
6. To renew your support license, click on the renew button and follow the instructions to extend your support license.

LICENSE INFO	
Office Name	Mississauga Smiles Dentistry
Email	agobin@smartreader.ca
License Number	XXXX-XXXX-XXXX-1111
Support Expiry Date	Saturday January 02, 2021

Support license expires in 344 day(s).

[Click here to renew support license](#)

LICENSE INFO	
Office Name	Mississauga Smiles Dentistry
Email	agobin@smartreader.ca
License Number	XXXX-XXXX-XXXX-1111
Support Expiry Date	Wednesday February 12, 2020

Support license expires in 19 day(s).


[Click here to renew support license](#)

3.4 Setup Autoclave Information

On this page you can add, edit or delete your autoclaves.

3.4.1 Add Autoclave

1. To add a new autoclave, please enter the preferred name for your autoclave. (Example Midmark A. Midmark 9, e.t.c) - This is the name your office prefers to call the autoclave.


AUTOCLAVE SETUP

Add new autoclave To add a new autoclave, please have your autoclave serial number ready.

Name Manufacturer

Serial Number Model

Select Type Log Directory

Modify Autoclave

Autoclave Name	Manufacturer	Model	Serial Number	Type	Log Directory
----------------	--------------	-------	---------------	------	---------------

2. Choose the Manufacturer and Model

Add new autoclave To add a new autoclave, please have your autoclave serial number ready.

Name Manufacturer

Serial Number Model

Select Type Log Directory

Modify Autoclave

Autoclave Name	Manufacturer	Model	Serial Number
----------------	--------------	-------	---------------

Select manufacturer

- Select manufacturer
- Scican
- Midmark
- Ritter
- Pelton-Crane
- W&H Lisa
- Flight Dental
- Masteri
- Tuttnauer
- Other

Select model

- Select model
- Select model
- M9
- M11
- M3
- M9 1st Gen
- M11 1st Gen
- Other

3. Enter the serial number of your autoclave (Please find it in the autoclave's manufacturer manual, manufacturer documentation or on the device)
4. Choose the type of autoclave (Regular or Class B)

Add new autoclave To add a new autoclave, please have your autoclave serial number ready.

Name Manufacturer

Serial Number Model

Select Type Log Directory

Modify Autoclave

Autoclave Name	Manufacturer	Model	Serial Number	Type	Log Directory
----------------	--------------	-------	---------------	------	---------------

5. If Class B chosen, the Bowie Dick test becomes activated in the load management.
6. The Log directory is the location of where the logs are stored. With our Smart Reader wireless office setup, the directory where the logs of an autoclave are stored can be linked to the autoclave in the software. If set, the software automatically reads the most recent log file whenever the autoclave is selected in the sterilization stage. If the location is known, click on the browse button to select the log directory. If your office doesn't have our Smartreader wireless office setup, leave this field black.
7. After all the input has been added and verified, click on the **Add Autoclave** button to add the autoclave.
8. You will receive a notification that the autoclave has been successfully added.

Autoclave added successfully

Add new autoclave To add a new autoclave, please have your autoclave serial number ready.

Name Manufacturer

Serial Number Model

Select Type Log Directory

Modify Autoclave

Autoclave Name	Manufacturer	Model	Serial Number	Type	Log Directory
Midmark A	Midmark	M9	129765377	Regular	C:\Users\Momin\Desktop\Agobin\Logs\Midmark

3.4.2 Missing Autoclave Manufacturer and/or Model

1. If your autoclave's manufacturer or model is not in the list you can add it to the list by choosing other in the manufacturer or model list. After selecting other, an input field would appear, enter the manufacture and/or model.

3.4.3 Edit Autoclave

If the information for your autoclave is not correct, you can edit the information

1. Select the autoclave from the list.
2. Enter the new information. The name, serial number, manufacturer, model, type and log directory can be changed.
3. Click on the update button to update the autoclave's information.

3.4.4 Delete Autoclave

If your autoclave is not in use anymore, you can delete the autoclave.

1. Select the autoclave from the list you want to delete.
2. Click on the delete button.
3. You will be prompted to confirm if you want to delete the autoclave.
4. Click yes to confirm.
5. The autoclave will be deleted and removed from the list.
6. You will receive a notification that autoclave has been successfully added.
7. Though the autoclave has been deleted, sterilization reports of the autoclave will still be available.

Add new autoclave

To add a new autoclave, please have your autoclave serial number ready.

Name Manufacturer

Serial Number Model

Select Type Log Directory

Modify Autoclave

Confirmation

Do you want to delete?

Autoclave Name	Manufacturer	Model	Type	Log Directory
Midmark A	Midmark	M9	Regular	C:\Users\Momin\Desktop\Agobin\Logs\Midmark
Scican A	Scican	Statim 7000	Regular	
Smartclave A	Smartclave	L100	Class B	

3.5 Setup Operator Information

On this page you will add, delete or change information for your operators.

3.5.1 Add Operator

1. To add operator, please enter first name and last name.
2. Enter the PIN for the operator or click on the generate PIN button to generate a unique PIN for the operator.
3. To make an operator an administrator (this will give the operator administrator privileges to change information in the admin panel), click on the “Is Operator Admin?” checkbox.

Add new operator

Please provide the operator his/her PIN.

First Name

Enter PIN Generate PIN

Last Name

☐ Is Operator Admin?

Add Operator

Modify Operator

Update Delete Cancel

First Name	Last Name	PIN	Administrator Status
<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>

4. Click the Add Operator button.
5. You will receive a notification that operator has been successfully added.
6. You have completed the operator entry.
7. Please add all operators that will be using the system.
8. If an operator is an administrator, please provide pin to the operator for the operator to be able to access the admin panel of the application.

3.5.2 Edit Operator

You can edit operators selecting the operator from the list.

1. Select the operator you wish to edit on the list.

Add new operator

Please provide the operator his/her PIN.

First Name	<input type="text" value="Mary"/>	Enter PIN	<input type="text" value="6945"/>	<input type="button" value="Generate PIN"/>
Last Name	<input type="text" value="Jane"/>	<input type="checkbox"/> Is Operator Admin?		
<input type="button" value="Add Operator"/>				

Modify Operator

First Name	Last Name	PIN	Administrator Status
John	Doe	1234	1
Mary	Jane	6945	0

2. Click on the field you wish to change.
3. Enter the new information.
4. Once you have completed all the changes, click on the update button.
5. You will receive a notification confirming that the operator has been updated.

User updated successfully

3.5.3 Delete Operator

To delete an operator, select the operator you wish to delete from the list.

1. Select the operator you wish to delete on the list.

Add new operator

Please provide the operator his/her PIN.

First Name Mary

Enter PIN 6945

Generate PIN

Last Name Jane

☐ Is Operator Admin?

Add Operator

Modify Operator

Update

Delete

Cancel

First Name	Last Name	PIN	Administrator Status
John	Doe	1234	1
Mary	Jane	6945	0

2. Click on the Delete button.
3. You will be asked for confirmation if you want to proceed to delete the operator.
4. Click yes to delete the operator.
5. The operator will be deleted and you will receive a confirmation that the operator was successfully deleted.

User deleted successfully.

6. Though an operator has been deleted, all the sterilization records linked to that particular operator will still be available in the reports section.

3.6 Manage Packages

3.6.1 Creating new packages.

To create a new package:

1. Select instrument names from the instruments list and click on add. If an instrument is not in the list, enter the instrument name in the other instrument name box and click on the **Add** button. Adding instruments to the package is optional.
2. To remove an instrument from the package instruments list, select the instrument name and click on the **Remove** button.
3. To remove all instruments from the package instruments list, click on the **Reset List** button.
4. Adding instruments is an optional procedure. Packages can be created with just the package name, not adding any instruments.
5. Enter the package name.
6. Click on **Save Package** button to save the package.
7. You will receive a notification that the package was successfully created.

☺ Create new package

Select Instruments to Add

Blades

Burs

Curettes

Carver

Cord Packers

Crown Removers

Crown Seaters

Drill bits

Dental Needles

Flawers

Add Other Instrument

Name

Add

Instruments in Package

Burnisher
Calipers
Cutters
Condenser

Add

Remove

Reset List

Package name Exam1

Save Package

3.6.2 Modifying existing packages.

To update a package:

1. Select the package name from the select package drop down box.
2. To add instruments to the package, select the instruments from the Instruments to Add list and click on the Add button.
3. If an instrument is not on the list, enter the instrument name in the other instrument name textbox and click on the Add button to add the instrument to the list.
4. To remove an instrument from the list, select the instrument(s) from the Instruments in selected package list and click on the Remove button to remove the package.
5. To remove all instruments from the package list, click on the Remove all instruments button.
6. After modifying the package, click on the update button.
7. You will receive a notification that the package has been updated.

To delete a package:

1. Select the package name from the select package drop down box.
2. Click on the Delete button.
3. You will be prompted for confirmation to delete the selected package.
4. Click yes, to continue.
5. You will receive a notification that the package has been successfully deleted.

☐ Create new package

☐ Modify Package

Select Package Bib Clips Delete Update

Select Instruments to Add

- Blades
- Calipers
- Cutters
- Carver
- Condenser
- Cord Packers
- Crown Removers
- Crown Seaters
- Drill bits
- Dental Needles

Add

Remove

Add Other Instrument

Name

Add

Instruments in Selected Package

- Burnisher
- Burs
- Curettes

Remove all instruments

3.6.3 Viewing packages.

To view a package:

1. Select the package name in the package list.
2. Instruments in the selected package, if any, would be displayed in the instruments in package list box.

☑ Create new package

☑ Modify Package

☑ View Packages

Select Package to View

Bib Clips
Burs
Cavitron Tips
Crown Bridge
Earbuds
Endo
Examination
Handpieces
Hygiene
Implant Scalers

Instruments in Bib Clips

Burnisher
Blades
Curettes

3.7 Backup and Restore

3.7.1 Creating a backup

To create a backup of the application data:

1. Click on the Backup Data button.
2. A dialog box would open, choose the location you want to store the backup file in.
3. The system generates a name for the backup file, change the name to something else if you want to.
4. Click on the save button to save the backup in the selected location.
5. You will receive a notification that the backup was successfully created.

3.7.2 Setting up automatic backup

To set up automatic backup:

1. Click on the browse button.
2. Select the folder where you wish to store the automatic backups.
3. Click on the Ok button to choose the selected folder.
4. Click on the Update Backup Location button to set the automatic backup.
5. You will receive a notification that automatic backup has been set.


To change the automatic backup location:

1. Click on the browse button.
2. Select the folder/directory you want to change to.
3. Click on the Ok button to choose the selected folder.
4. Click on the Update Backup Location button to set the new backup location.
5. You will receive a notification that automatic backup has been set.

3.7.3 Restoring data

To restore application data:


1. Click on the Restore Data button.
2. Data restoration deletes all the application data and replaces it with the one in the backup file. This process is not reversible unless you created a backup before restoring the data. You will be prompted to provide confirmation if you want to proceed.
3. Click yes to proceed.
4. Select the backup file you wish to restore data from.
5. When the file is selected, click on the Open button.
6. You will receive a notification that data is being restored.
7. Upon successfully completion, you will receive another notification that data was successful restored.



3.8 Audit


In this section, you will be able to generate audit reports within a specified date range for your packages, autoclaves, operators (sterilization and biological indicator test operators) and reports based on various test (cycle temperature, cycle pressure, chemical indicator, biological indicator, type one, type four and type five tests) pass or fail status. To generate an audit report:

1. Click on the Audit menu on the navigation bar.
 2. Select the start date and the end date.
- A. To generate audit report for a package within the specified date range:
1. Select the package from the list.
 2. Click on the Get Report button.




The screenshot shows the 'AUDIT' section of the Sterilization Portal. At the top, there is a header bar with a bar chart icon and the word 'AUDIT'. Below this, there are two date selection fields: 'Select Start Date' with a calendar icon showing 'January 8, 2020' and '15', and 'Select End Date' with a calendar icon showing 'January 15, 2020' and '15'. Underneath the dates, there are five tabs: 'Packages', 'Autoclave', 'Sterilization Operator', 'Biological Indicator Operator', and 'Test Status'. The 'Packages' tab is currently selected. Below the tabs, there is a 'Select package' dropdown menu and a 'Get Report' button.

- B. To generate audit report for an autoclave within the specified date range:
1. Select the autoclave from the list.
 2. Click on the Get Report button.



The screenshot shows the 'AUDIT' section of the Sterilization Portal. At the top, there is a header bar with a bar chart icon and the word 'AUDIT'. Below this, there are two date selection fields: 'Select Start Date' with a calendar icon showing 'January 8, 2020' and '15', and 'Select End Date' with a calendar icon showing 'January 15, 2020' and '15'. Underneath the dates, there are five tabs: 'Packages', 'Autoclave', 'Sterilization Operator', 'Biological Indicator Operator', and 'Test Status'. The 'Autoclave' tab is currently selected. Below the tabs, there is a 'Select Autoclave' dropdown menu and a 'Get Report' button.

- C. To generate audit report for a sterilization operator within the specified date range:
1. Select the operator from the list.
 2. Click on the Get Report button.

 **AUDIT**

Select Start Date Select End Date

Packages

Autoclave


Sterilization Operator

Biological Indicator Operator

Test Status

Select Sterilization Operator

- D. To generate audit report for a biological indicator test operator within the specified date range:
1. Select the operator from the list.
 2. Click on the Get Report button.

 **AUDIT**

Select Start Date Select End Date

Packages

Autoclave

Sterilization Operator



Biological Indicator Operator

Test Status

Select B. I. Test Operator


- E. To generate audit report based on the various tests pass or fail status:
1. Select the test status, either **pass** or **fail**.
 2. Click on the Get Report button.

AUDIT

Select Start Date January 8, 2020  Select End Date January 15, 2020 

Packages Autoclave Sterilization Operator Biological Indicator Operator **Test Status**

Select B. I. Test, Cycle Temperature, Cycle Pressure Status.

☒ **Pass** ☐ **Fail** 

After the audit report has been generated, to sign the report enter the manager's initials in the appropriate text field. To email the report, enter the email(s) in the email text field and click on the Email Report button. To save the generated report locally, for example in a flash drive, click on the Save Report button. A save file dialog box will open, navigate to the location you want to store the report, enter the name of the file and click on the Save button. The report will be saved in the selected location.

STERILIZATION PORTAL
Mississauga smiles dentistry

AUDIT

Select Start Date: January 8, 2020 | Select End Date: January 15, 2020

Packages | Autoclave | Sterilization Operator | Biological Indicator Operator | Test Status

Select package: Bib Clips | Get Report

Number of packages sterilized: 1
Total number of cycles: 1

Enter manager initials to sign report:

Enter email addresses: john.doe@smartdental.ca

To email report to multiple emails, separate emails with a comma. E.g abc@gmail.com, def@hotmail.com

Email Report | Save Report

Click the Save Report button to save report

Click the Email Report button to email report

Package Information			
Bib Clips(Unwrapped)(1)			
Clean Method	Hydrex	Cycle No	7105
Package ID(s)	111121212121	Cycle Temperature	Fail-134C
Autoclave Name	Midmark A	Cycle Pressure	Pass-2026PS
Autoclave Manufacturer	Midmark	Sterilization Duration	3 Minutes
Autoclave Model	M9	Chemical Test/Lot#	Pass-2341
Autoclave SN	7823723	Cycle Time	15:03
Bowie Dick Test	N/A		
Cycle Date	01/14/2020		
Cycle Operator	John Doe		

Biological Indicator Test Information

B. I. Spot Test Result	Pass	B. I. Test Lot#	1231
Type 1 Test Status	Pass	Type 4 Test Status	Pass
Type 2 Test Status	Pass	Type 5 Test Status	Pass

You have completed the Administration setup. Please log out of the Administration setup by clicking on the Logout button on the left side of the screen. You will be returned to the home page.

The Administrator can re-log in to the administration page for changes or updates by clicking on the Administration button on the left navigation pane.

3.9 Logging Out from Administration Panel


To logout from the admin panel:

1. Click on the Logout button.
2. You will be redirected to the Sterilization page.

4 Main Menu

4.1 Login

To log into the application the operator will need the office ID and PIN. That information was sent to the administrator.




Office Login

Office ID

PIN

Login

[Forgot PIN?](#)



STERILIZATION PORTAL
Mississauga Smiles Dentistry

[SIGN OUT](#) **STERI SMART** Wed, 15 Jan 2020 11:30

STERILIZATION

Select Cleaning Method

Select Package, Number and Type Present
Enter/Scan Package Barcode(s) (separate barcodes with commas if multiple)

Number	Type	Number	Type	Number	Type	Number	Type	Number	Type	Number	Type
Bib Clips		Burs		Cavitron Tips		Crown Bridge		Earbuds		Endo	

Sterilization

Select Operator

Select Autoclave

Cycle Number

Cycle Date

Sterilization Time

Enter values for temperature and pressure with their units. C/F for temperature and PSI/KPA for pressure.

☐ Maximum Temperature ☐ Maximum Pressure

☐ Chemical Indicator Pass

4.2 Navigation

You can use the links on the left side of the screen for navigation.

Please note that you can only see the report without the B. I. test, but to complete the report you will need the B. I. test for the autoclave completed.

No complete reports will be compiled and saved without the B. I. Test

4.2.1 Step 1 – Printing of Labels

In this step, the application generates and prints barcode labels which will be applied to all the packages before they go into the autoclave.

To print barcode labels:

1. Click on the print labels option in the navigation pane.
2. Click on the print labels button.
3. Make sure the right printer is selected.
4. Enter the number of copies to print.
5. Uncheck the “Fit picture to frame” checkbox.
6. Finally, click on the print button.

The desired number of barcode labels would be printed. Apply the printed barcode labels to the packages.



PRINT LABELS

Print Labels

Click the print button to print labels

Print Labels



4.2.2 Step 2 – Sterilization

In this step, you report how the instruments were cleaned (options include dental washer, manual, ultrasonic and hydrim). You also report the type of packages (load), the package ID(s), the number of each package present and the type of each package. Next, you input all information on the sterilization cycle (operator, autoclave, cycle number, temperature, pressure, chemical indicator test and lot number)

To enter the sterilization cycle information:

1. Click on the sterilization link on the left navigation pane if the sterilization page isn't already opened.
2. Select the cleaning method by clicking on the dropdown menu and select the method which applied to how the instruments were cleaned.
3. After selecting the cleaning method, enter or scan the package barcode labels. If all the packages have the same label, scan/enter the package ID only once. If you are using pre-printed barcode labels (or if some packages have different barcode/package IDs) scan each different barcode label separating them with a comma (",").
4. Next, select all packages present. Select the number present and the package type.
5. Select the sterilization operator from the dropdown list.
6. Select the Autoclave used in sterilization. If the Log directory of the autoclave has been set and the directory exists, the system would read the most recent log file from the directory and capture the needed information and fill in the required fields.
7. If the log directory of the autoclave wasn't set and the location of the cycle log file is known (if the cycle log file is in a USB or SD card, insert the USB/SD card before you proceed), click the "Upload Log" button to upload the log file. After clicking the button, a dialog opens, navigate to the location where the cycle log file is found. Select the log file and click the open button. The system reads the log file, extracts and fill all the needed information.
8. If the log directory of the autoclave wasn't set and the cycle log file doesn't exist, manually enter the cycle number, cycle date, hour, minute, sterilization time (optional field).
 - a. If the maximum temperature for sterilization was reached, check the maximum temperature checkbox and enter the temperature if known (entering the maximum temperature value is optional).
 - b. If the maximum pressure for sterilization was reached, check the maximum pressure checkbox and enter the maximum pressure (entering the maximum pressure value is optional).
 - c. If the chemical indicator test is pass, check the chemical indicator pass checkbox and enter the chemical indicator lot number.

9. Verify that all input data is correct.
10. Click the submit button to finally submit all information.
11. You will receive a notification that data was submitted successfully.

Data saved successfully

4.2.2.1 Smart Reader Data Logger

If you have a Smart Reader Data Logger, it is easy to upload the data.

1. Once the cycle is complete, please remove the SD card from the Smart Reader
2. Insert the SD card into the PC
3. Open a new sterilization record by clicking the sterilization link in the navigation pane.
4. Click on the Upload Log button.

STERILIZATION

Select Cleaning Method

Select Package, Number and Type Present

Enter/Scan Package Barcode(s) (separate barcodes with commas if multiple)

Bib Clips: Number Type
 Burs: Number Type
 Cavitron Tips: Number Type
 Crown Bridge: Number Type
 Earbuds: Number Type
 Endo: Number Type

Sterilization

Select Operator

Select Autoclave

Cycle Number


Cycle Date Hour Minute

Sterilization Time

Enter values for temperature and pressure with their units. C/F for temperature and PSI/KPA for pressure.

☐ Maximum Temperature ☐ Maximum Pressure
☐ Chemical Indicator Pass

- File explorer will open and choose your cycle log file.



- Double click the cycle
- The data will populate the fields. (Cycle number, temperature, pressure, sterilization duration, date and time)
- Continue with the chemical indicator test data entry

4.2.3 Step 3 – Biological Indicator Test

In this step you input all information on the biological indicator test. To enter the cycle biological indicator test information:

1. Click on the B. I. Test link in the left navigation pane.
2. Select the cycle with pending biological indicator test result.
3. Select the biological indicator test operator from the drop down list.
4. Enter the biological indicator lot number.
5. Check the biological indicator type 1, 4, 5 checkboxes if all test pass. By default, the type 1, 4 and 5 chemical test are all set to pass (green in color). If any of the test fail, check the respective checkbox, it changes color to red.
6. Verify that all data input are correct.
7. Click the submit button to submit the biological indicator test result.
8. You will receive a notification that the biological indicator test result was submitted successfully.

STERILIZATION PORTAL
Mississauga Smiles Dentistry

BIOLOGICAL INDICATOR TEST

Biological Indicator Test Result

Select B. I. Test Operator: Jane Doe

Select Autoclave: Midmark A

Select B. I. Test Date: August 1, 2019

B. I. Spot Test Result: Pass

Enter B. I. Lot Number: 4321

☒ Type 1 Test Pass ☒ Type 4 Test Pass ☒ Type 5 Test Pass

Submit

Cycle Operator	Cycle Number	Autoclave	Manufacturer	Cycle Date	BI Result
John Doe	7106	Midmark A	Midmark	August 1, 2019	PENDING

Biological Test result saved successfully

4.2.4 Step 4 – Reports

Sterisart creates sterilization reports that are compliant with health regulations and it provides traceability report of sterilization cycle records. The Sterilization Portal will store your results and attach them to the cycle report for the autoclave. You can generate cycle reports daily, by date range, by autoclave, by sterilization user, by biological indicator test user and by package ID (barcode). All reports are formatted to simplify the overview and have names and signatures of operators and managers.

To generate a cycle report:

1. Click on the reporting link in the left navigation pane.
2. To generate daily reports:
 - a. Click the daily tab under sterilization report.
 - b. Select the date you wish to generate the report.
 - c. Click the generate daily report button.
 - d. Report would be generated if any cycles exists on that particular date.



REPORTING

Sterilization Report

Daily

Date Range

Autoclave

Sterilization User

Biological Indicator Test User

Barcode

Generate Daily Reports

Select Date

January 14, 2020



Generate Daily Report

3. To generate reports by date range:
 - a. Click the date range tab under the sterilization report.
 - b. Select the start (from) date.
 - c. Select the end (to) date.
 - d. Click the generate date range report button.
 - e. Reports would be generated if any cycles exists within the selected dates.



REPORTING

Sterilization Report

Daily

Date Range

Autoclave

Sterilization User

Biological Indicator Test User

Barcode

Generate Date Range Reports

From

January 1, 2020




To

January 15, 2020



Generate Date Range Report

4. To generate reports by autoclave:
 - a. Click on the autoclave tab under sterilization report.
 - b. Select the autoclave from the dropdown list.
 - c. To view all the cycle reports of the selected autoclave, click on the generate all report button.
 - d. To view all the cycle reports of the selected autoclave on a particular date, select the desired **from** date and click on the generate date report button.
 - e. To view all the cycle reports of the selected autoclave by date range, select the **start** (from) date and the **end** (to) date and click on the generate date report button.

 **REPORTING**

Sterilization Report

Daily

Date Range

Autoclave

Sterilization User

Biological Indicator Test User

Barcode

Generate Reports by Autoclave

Select Autoclave

Midmark A

Generate All Report

Filter Autoclave Report by Date

From

January 15, 2020

15

To (Optional)

January 15, 2020

15

Generate Date Report

5. To generate reports by sterilization user:
 - a. Click on the sterilization user tab under sterilization report.
 - b. Select the user from the dropdown list.
 - c. To view all reports done by the selected sterilization user, click on the generate all report button.
 - d. To view all cycle reports of the selected sterilization user done on a particular date, select the date (from) and click on the generate date report button.
 - e. To view all the cycle reports of the selected sterilization user by date range, select the start (from) date and the end (to) date and click on the generate date report button.



REPORTING

Sterilization Report

Daily

Date Range

Autoclave

Sterilization User

Biological Indicator Test User

Barcode

Generate Reports by Sterilization User

Select User Jane Doe

Generate All Report

Filter Sterilization User Report by Date

From January 15, 2020

15

To (Optional) January 15, 2020

15

Generate Date Report

6. To generate reports by biological indicator test user:
 - a. Click on the biological indicator test user tab under sterilization report.
 - b. Select the user from the dropdown list.
 - c. To view all reports done by the selected biological indicator test user, click on the generate all report button.
 - d. To view all cycle reports of the selected biological indicator test user done on a particular date, select the date (from) and click on the generate date report button.
 - e. To view all the cycle reports of the selected biological indicator test user by date range, select the start (from) date and the end (to) date and click on the generate date report button.



REPORTING

Sterilization Report

Daily

Date Range

Autoclave

Sterilization User

Biological Indicator Test User

Barcode

Generate Reports by Biological Indicator Test User

Select User

John Doe

Generate All Report

Filter Biological Indicator User Report by Date

From

January 15, 2020

15

To (Optional)

January 15, 2020

15

Generate Date Report

7. To generate report by barcode:
 - a. Click on the barcode tab under sterilization report.
 - b. Enter or scan the barcode ID number.
 - c. Click the generate report button.



REPORTING

Sterilization Report

Daily

Date Range

Autoclave

Sterilization User

Biological Indicator Test User

Barcode

Generate Reports by Barcode

Enter or scan barcode

Generate Report

Generate cycle reports can be stored as a PDF file. To save the generated report:

1. After the sterilization report has been generated, to save the report, click on the save report button.
2. A dialog opens, navigate to the directory where you wish to store the PDF report.
3. Enter the name of the report if you wish to change the default name.
4. After verifying the directory and the name of the PDF file, click on the save button.
5. A PDF report is generated and opened by the default PDF file viewer in the system.
6. To email the report, enter the email(s) in the email box, if multiple emails, separate emails with a comma. After entering all emails, click on the email report button. You will receive a confirmation message that email has been sent.

Note: All sterilization cycles with a pending biological indicator test are excluded from the results. Meaning only cycles with either a pass or fail biological indicator test are saved in the PDF report.


To sign a report as the manager:

1. After generating the report, enter the manager's initials in the box that requests for the manager's initials.
2. Click the save report button and continue as outlined above to save or email the PDF report.

4.2.5 Administration

Administration of the portal. You need to have admin rights to access this option. To change system configuration (i.e to change office information, autoclave, operators, packages and/or backup settings):

- a. Click on the Administration link in the navigation pane.
- b. A dialog opens requesting you to enter your PIN number.
- c. Enter your PIN number and click continue.
- d. If admin PIN validation is successful, the user would be directed to the admin portal to add, remove or modify system configuration.



5 Help Center

All information on how to get help, documentation and manuals.

5.1 Help

The telephone for support is (905) 997-7171. The email is support@smartsoftware.ca

5.2 Training

For additional training please contact the office at (905) 997-7171

5.3 FAQ

Here will be all user questions and answers by our technical department in regard to most common questions.

5.3.1 Documentation

All Sterilization information like IPAC info, RCDSO info and Public health information will be available here. MSDS sheets will be added as the database receives email from users.

5.3.2 Manual

The PDF Manual is located in this section.

6 Sign Out

To sign out from the system:

1. Click the red sign out button located at the top right hand side corner of the screen.
2. You will be redirected to the login screen.

7 Appendix

7.1 LICENSE AND SERIAL NUMBER

To use our Sterilization Portal, please use the following information:

Serial Number:

License Number:

Please use this information to register your office in the system.